


Ozsmart is an active therapeutic medical device designed for ozone therapy applications in clinical and outpatient settings, providing an effective solution for the requirements of **modern ozone therapy**.

The biochemical action of ozone in the body supports blood circulation in all tissues, including the central nervous system, exhibits bactericidal, fungicidal, and virucidal properties, has an immunorestorative effect, contributes to cellular energy processes, and promotes tissue regeneration.



- Mobile and ergonomically designed compact solution
- 7" touchscreen display
- Integrated RFID card reader for secure user authentication
- Advanced nozzle flow control system
- Multilingual user interface with operator guidance
- Catalytic ozone destructor
- Comprehensive automatic self-test at every device start-up



Ozsmart offers the full spectrum of ozone therapy procedures:

- Major autohemotherapy (MAH)
- Minor autohemotherapy (MiAH)
- Gas collection for rectal/vaginal insufflation
- Filling of treatment bags
- Ozonation of liquids

Indications

- Inflammatory conditions (bacterial, viral, and fungal)
- Ischemic conditions
- Weakened immune system
- Chronic pain

Technical Parameters

Device Name	Ozsmart
Power Supply	230 V / 50 Hz
Power Consumption	100 VA
Ozone Concentration	1 μ g / Nml - 80 μ g
Flow Rate	0,2 or 0,8 l/min. depending on mode
Dimensions (W x H x D)	380 x 310 x 135 mm
Dimensions of the Control and Safety Interlock Unit	150x110x190 mm
Weight of the Main Unit	13,6 kg
Product Category	II b

Manufacturer:

ZAT a.s.
K Podlesí 541
26101 Příbram
Czech Republic

Distributor:

MEDIATRADE s.r.o.
Polična No.34
75701 Polična
Czech Republic

Contact person:

Aleksandr Rudometov
Phone: +420 736 780 608
Email: rudometov@mediatrade.cz
www.mediatrade.cz

